


CRISTALIZOARE

CRISTALIZOR

- Este confectionat din stila DURAN
- Poate fi autoclavat



Cristalizer cu cioc

Cod	Produs	Volum	Diametrul exterior	Inaltime	Mod de ambalare
CRTC3975.0010	Cristalizer cu cioc	20 ml	40 mm	25 mm	10 bucati
CRBC3976.0010	Cristalizer cu cioc	40 ml	50 mm	30 mm	10 bucati
CRFC3977.0010	Cristalizer cu cioc	60 ml	60 mm	35 mm	10 bucati
CRSC3978.0010	Cristalizer cu cioc	100 ml	70 mm	40 mm	10 bucati
CRQC3979.0010	Cristalizer cu cioc	150 ml	80 mm	45 mm	10 bucati
CRPC3980.0010	Cristalizer cu cioc	300 ml	95 mm	55 mm	10 bucati
CRZC3981.0001	Cristalizer cu cioc	500 ml	115 mm	65 mm	bucata
CRGC3982.0001	Cristalizer cu cioc	900 ml	140 mm	75 mm	bucata

CRFC3983.0001	Cristalizer cu cioc	2000 ml	190 mm	90 mm	bucata
CRYC3984.0001	Cristalizer cu cioc	3500 ml	230 mm	100 mm	bucata

Cristalizer fara cioc

Cod	Produs	Volum	Diametrul exterior	Inaltime	Mod de ambalare
CRMC3985.0010	Cristalizer fara cioc	20 ml	40 mm	25 mm	10 bucati
CRIC3986.0010	Cristalizer fara cioc	40 ml	50 mm	30 mm	10 bucati
CRTC3987.0010	Cristalizer fara cioc	60 ml	60 mm	35 mm	10 bucati
CRRC3988.0010	Cristalizer fara cioc	100 ml	70 mm	40 mm	10 bucati
CRDC3989.0010	Cristalizer fara cioc	150 ml	80 mm	45 mm	10 bucati
CRHC3990.0010	Cristalizer fara cioc	300 ml	95 mm	55 mm	10 bucati
CRJC3991.0001	Cristalizer fara cioc	500 ml	115 mm	65 mm	bucata
CRSC3992.0001	Cristalizer fara cioc	900 ml	140 mm	75 mm	bucata
CRXC3993.0001	Cristalizer fara cioc	2000 ml	190 mm	90 mm	bucata
CRNC3994.0001	Cristalizer fara cioc	3500 ml	230 mm	100 mm	bucata