

BALOANE COTATE



BALON COTAT CLASA A

- Dopul poate fi din sticla sau din plastic
- Este confectionat din sticla borosilicata 3,3 mm, termorezistenta, transparenta sau bruna

Balon cotat clasa A, transparent, cu dop din plastic

Cod	Produs	Toleranta (± ml)	SN	Mod de ambalare
BRTC4406.0002	Balon cotat, clasa A, 5 ml	0,025	7/16	2 bucati
BRGC4407.0002	Balon cotat, clasa A, 5 ml	0,025	10/19*	2 bucati
BRHC4408.0002	Balon cotat, clasa A, 10 ml	0,025	7/16	2 bucati
BRBC4409.0002	Balon cotat, clasa A, 10 ml	0,025	10/19*	2 bucati
BRMC4410.0002	Balon cotat, clasa A, 20 ml	0,040	10/19	2 bucati
BRQC4411.0002	Balon cotat, clasa A, 25 ml	0,040	10/19	2 bucati
BRZC4412.0002	Balon cotat, clasa A, 50 ml	0,060	12/21	2 bucati
BRLC4413.0002	Balon cotat, clasa A, 50 ml	0,060	14/23*	2 bucati
BRPC4414.0002	Balon cotat, clasa A, 100 ml	0,100	12/21	2 bucati
BROC4415.0002	Balon cotat, clasa A, 100 ml	0,100	14/23*	2 bucati

BRDC4416.0002	Balon cotat, clasa A, 200 ml	0,150	14/23	2 bucati
BRSC4417.0002	Balon cotat, clasa A, 250 ml	0,150	14/23	2 bucati
BRCC4418.0002	Balon cotat, clasa A, 500 ml	0,250	19/26	2 bucati
BRZC4419.0002	Balon cotat, clasa A, 1000 ml	0,400	24/29	2 bucati
BRNC4420.0002	Balon cotat, clasa A, 2000 ml	0,600	39/32	2 bucati
BRRC4421.0001	Balon cotat, clasa A, 5000 ml	1,200	34/35	bucata

*

NOTA: = gat larg

Balon cotat clasa A, transparent, cu dop din sticla

Cod	Produs	Caracteristici slif (Ø x h) mm	Mod de ambalare
BMPE116.0001	Balon cotat, clasa A, 5 ml	7/11	bucata
BMIE117.0001	Balon cotat, clasa A, 10 ml	7/11	bucata
BMBE118.0001	Balon cotat, clasa A, 25 ml	10/13	bucata
BMSE119.0001	Balon cotat, clasa A, 50 ml	10/13	bucata
BMDE120.0001	Balon cotat, clasa A, 100 ml	12/14	bucata
BMLE121.0001	Balon cotat, clasa A, 200 ml	14/15	bucata
BMOE122.0001	Balon cotat, clasa A, 250 ml	14/15	bucata
BMAE123.0001	Balon cotat, clasa A, 500 ml	16/16	bucata
BMZE124.0001	Balon cotat, clasa A, 1000 ml	19/17	bucata
BMHE125.0001	Balon cotat, clasa A, 2000 ml	24/20	bucata

Balon cotat clasa A, brun, cu dop din plastic

Cod	Produs	Toleranta (± ml)	SN	Mod de ambalare
BRYC4422.0002	Balon cotat, clasa A, 5 ml	0,025	7/16	2 bucati
BRXC4423.0002	Balon cotat, clasa A, 5 ml	0,025	10/19*	2 bucati
BRUC4424.0002	Balon cotat, clasa A, 10 ml	0,025	10/19*	2 bucati
BRJC4425.0002	Balon cotat, clasa A, 20 ml	0,040	10/19	2 bucati
BRHC4426.0002	Balon cotat, clasa A, 25 ml	0,040	10/19	2 bucati
BRFC4427.0002	Balon cotat, clasa A, 50 ml	0,060	12/21	2 bucati
BRTC4428.0002	Balon cotat, clasa A, 50 ml	0,060	14/23*	2 bucati
BRBC4429.0002	Balon cotat, clasa A, 100 ml	0,100	12/21	2 bucati

BRHC4430.0002	Balon cotat, clasa A, 100 ml	0,100	14/23*	2 bucati
BRLC4431.0002	Balon cotat, clasa A, 200 ml	0,150	14/23	2 bucati
BRMC4432.0002	Balon cotat, clasa A, 250 ml	0,150	14/23	2 bucati
BRPC4433.0002	Balon cotat, clasa A, 500 ml	0,250	19/26	2 bucati
BRBC4434.0002	Balon cotat, clasa A, 1000 ml	0,400	24/29	2 bucati

*

NOTA: = gat larg

Balon cotat clasa A, brun, cu dop din sticla

Cod	Produs	Caracteristici slif (Ø x h) mm	Mod de ambalare
BMCE136.0001	Balon cotat, clasa A, 5 ml	7/11	bucata
BMZE137.0001	Balon cotat, clasa A, 10 ml	7/11	bucata
BMPE138.0001	Balon cotat, clasa A, 25 ml	10/13	bucata
BMLE139.0001	Balon cotat, clasa A, 50 ml	10/13	bucata
BMYE140.0001	Balon cotat, clasa A, 100 ml	12/14	bucata
BMEE141.0001	Balon cotat, clasa A, 200 ml	14/15	bucata
BMWE142.0001	Balon cotat, clasa A, 250 ml	14/15	bucata
BMFE143.0001	Balon cotat, clasa A, 500 ml	16/16	bucata
BMIE144.0001	Balon cotat, clasa A, 1000 ml	19/17	bucata
BMOE145.0001	Balon cotat, clasa A, 2000 ml	24/20	bucata

BALON COTAT CLASA B

- Este confectionat din sticla borosilicata 3,3 mm, termorezistenta, transparenta
- Dopul este confectionat din sticla

Cod	Produs	Caracteristic slif (Ø x h) mm	Mod de ambalare
BMQE156.0001	Balon cotat, clasa B, 5 ml	7/11	bucata
BMZE157.0001	Balon cotat, clasa B, 10 ml	7/11	bucata
BMNE158.0001	Balon cotat, clasa B, 25 ml	10/13	bucata
BMFE159.0001	Balon cotat, clasa B, 50 ml	10/13	bucata
BMTE160.0001	Balon cotat, clasa B, 100 ml	12/14	bucata
BMPE161.0001	Balon cotat, clasa B, 200 ml	14/15	bucata

BMLE162.0001	Balon cotat, clasa B, 250 ml	14/15	bucata
BMCE163.0001	Balon cotat, clasa B, 500 ml	16/16	bucata
BMNE164.0001	Balon cotat, clasa B, 1000 ml	19/17	bucata
BMME165.0001	Balon cotat, clasa B, 2000 ml	24/20	bucata

BALON COTAT DIN PLASTIC

- Este confectionat din PMP (polimetilpentena)
- Dopul este din PP (polipropilena)
- Poate fi autoclavat
- Clasa B



Balon cotat din TPX

- rezistent (pentru scurt timp) la temperaturi de pana la + 175 °C

Cod	Produs	Capacitate	SN slif (Ø x h) mm	Inaltime	Mod de ambalare
BRTC3052.0001	Balon cotat, din plastic	25 ml	10/19	115 mm	bucata
BREC3053.0001	Balon cotat, din plastic	50 ml	10/19	150 mm	bucata
BRQC3054.0001	Balon cotat, din plastic	100 ml	14/23	180 mm	bucata
BRPC3055.0001	Balon cotat, din plastic	250 ml	19/26	235 mm	bucata

BRAC3056.0001	Balon cotat, din plastic	500 ml	19/26	270 mm	bucata
BRZC3057.0001	Balon cotat, din plastic	1000 ml	24/29	310 mm	bucata

Balon cotat din PP (polipropilena)

- rezistent (pentru scurt timp) la temperaturi de pana la + 140 °C

Cod	Produs	Capacitate	SN slif (Ø x h) mm	Inaltime	Mod de ambalare
BRLC3058.0001	Balon cotat, din plastic	25 ml	10/19	115 mm	bucata
BRFC3059.0001	Balon cotat, din plastic	50 ml	10/19	150 mm	bucata
BRJC3060.0001	Balon cotat, din plastic	100 ml	14/23	180 mm	bucata
BRNC3061.0001	Balon cotat, din plastic	250 ml	19/26	235 mm	bucata
BRGC3062.0001	Balon cotat, din plastic	500 ml	19/26	270 mm	bucata
BRCC3063.0001	Balon cotat, din plastic	1000 ml	24/29	310 mm	bucata