

BALOANE CU FUND PLAT / ROTUND

BALON CU FUND PLAT

- Confectionate din sticla borosilicata 3,3 termorezistenta
- SN= diametrul partii de imbinare/ inaltimea partii de imbinare (mm)



Cod	Produs	SN	Ø balon (mm)	H _{totala} (mm)	Mod de ambalare
BMVE476.0001	Balon fund plat 50 ml	14/23	53	88	bucata
BMNE477.0001	Balon fund plat 50 ml	19/26	53	94	bucata
BMJE478.0001	Balon fund plat 100 ml	14/23	65	100	bucata
BMIE479.0001	Balon fund plat 100 ml	19/26	65	106	bucata
BMFE480.0001	Balon fund plat 100 ml	24/29	65	109	bucata
BMTE481.0001	Balon fund plat 250 ml	19/26	88	134	bucata
BMCE482.0001	Balon fund plat 250 ml	24/29	88	139	bucata
BMYE439.0001	Balon fund plat 250 ml	29/32	88	145	bucata
BMAE483.0001	Balon fund plat 500 ml	19/26	110	158	bucata

BMBE484.0001	Balon fund plat 500 ml	24/29	110	161	bucata
BMJE440.0001	Balon fund plat 500 ml	29/32	110	167	bucata
BMSE486.0001	Balon fund plat 1000 ml	24/29	140	196	bucata
BMKE441.0001	Balon fund plat 1000 ml	29/32	140	199	bucata
BMXE442.0001	Balon fund plat 2000 ml	29/32	175	234	bucata

BALON CU FUND ROTUND

- Confectionat din sticla borosilicata 3,3 termorezistenta
- SN= diametrul partii de imbinare/ inaltimea partii de imbinare (mm)



Cod	Produs	SN	Ø balon (mm)	H _{totala} (mm)	Mod de ambalare
BMSE511.0001	Balon fund rotund 5 ml	10/19	27	56	bucata

BMVE370.0001	Balon fund rotund 5 ml	14/23	27	60	bucata
BMTE371.0001	Balon fund rotund 10 ml	10/19	33	66	bucata
BMBE372.0001	Balon fund rotund 10 ml	14/23	33	68	bucata
BMDE373.0001	Balon fund rotund 25 ml	10/19	41	70	bucata
BMIE374.0001	Balon fund rotund 25 ml	14/23	41	75	bucata
BMOE375.0001	Balon fund rotund 25 ml	19/26	41	79	bucata
BMZE376.0001	Balon fund rotund 50 ml	14/23	53	88	bucata
BMVE377.0001	Balon fund rotund 50 ml	19/26	53	94	bucata
BMPE378.0001	Balon fund rotund 50 ml	24/29	53	97	bucata
BMBE379.0001	Balon fund rotund 50 ml	29/32	53	105	bucata
BMWE380.0001	Balon fund rotund 100 ml	14/23	65	100	bucata
BMFE381.0001	Balon fund rotund 100 ml	19/26	65	106	bucata
BMME382.0001	Balon fund rotund 100 ml	24/29	65	109	bucata
BMNE383.0001	Balon fund rotund 100 ml	29/32	65	117	bucata
BMXE384.0001	Balon fund rotund 250 ml	19/26	88	134	bucata
BMEE385.0001	Balon fund rotund 250 ml	24/29	88	139	bucata
BMDE386.0001	Balon fund rotund 250 ml	29/32	88	145	bucata
BMRE387.0001	Balon fund rotund 500 ml	19/26	110	158	bucata
BMSE388.0001	Balon fund rotund 500 ml	24/29	110	161	bucata
BMJE475.0001	Balon fund rotund 500 ml	29/32	110	167	bucata
BMPE391.0001	Balon fund rotund 1000 ml	19/26	140	188	bucata
BMLE392.0001	Balon fund rotund 1000 ml	24/29	140	196	bucata
BMCE393.0001	Balon fund rotund 1000 ml	29/32	140	199	bucata
BMIE395.0001	Balon fund rotund 2000 ml	24/29	175	231	bucata
BMNE396.0001	Balon fund rotund 2000 ml	29/32	175	234	bucata
BMAE398.0001	Balon fund rotund 3000 ml	29/32	199	281	bucata
BMQE400.0001	Balon fund rotund 3000 ml	40/38	199	287	bucata
BMGE401.0001	Balon fund rotund 5000 ml	40/38	235	333	bucata
BMKE402.0001	Balon fund rotund 5000 ml	50/42	235	337	bucata
BMFE403.0001	Balon fund rotund 10000 ml	50/42	297	410	bucata
BMSE404.0001	Balon fund rotund 10000 ml	60/46	297	414	bucata

BALON CU FUND ROTUND, CU 2 GATURI

- Confectionat din sticla borosilicata 3,3 termorezistenta
- SN= diametrul partii de imbinare/ inaltimea partii de imbinare (mm)



Cod	Produs	SN-ul gatului central	SN-ul gatului lateral	H _{totala} (mm)	Mod de ambalare
BMNE512.0001	Balon fund rotund, 2 gaturi, 100 ml	19/26	14/23	106	bucata
BMME513.0001	Balon fund rotund, 2 gaturi, 100 ml	24/29	14/23	109	bucata
BMLE514.0001	Balon fund rotund, 2 gaturi, 250 ml	19/26	14/23	134	bucata
BMVE515.0001	Balon fund rotund, 2 gaturi, 250 ml	24/29	14/23	139	bucata
BMSE516.0001	Balon fund rotund, 2 gaturi, 250 ml	24/29	19/26	139	bucata
BMPE517.0001	Balon fund rotund, 2 gaturi, 250 ml	29/32	14/23	145	bucata
BMOE518.0001	Balon fund rotund, 2 gaturi, 250 ml	29/32	19/26	145	bucata
BMHE519.0001	Balon fund rotund, 2 gaturi, 500 ml	19/26	14/23	158	bucata
BMJE520.0001	Balon fund rotund, 2 gaturi, 500 ml	24/29	14/23	161	bucata
BMKE521.0001	Balon fund rotund, 2 gaturi, 500 ml	24/29	19/26	161	bucata
BMQE522.0001	Balon fund rotund, 2 gaturi, 500 ml	29/32	14/23	167	bucata
BMWE523.0001	Balon fund rotund, 2 gaturi, 500 ml	29/32	19/26	167	bucata

BMEE524.0001	Balon fund rotund, 2 gaturi, 1000 ml	24/29	14/23	196	bucata
BMTE525.0001	Balon fund rotund, 2 gaturi, 1000 ml	24/29	19/26	196	bucata
BMYE526.0001	Balon fund rotund, 2 gaturi, 1000 ml	24/29	24/29	196	bucata
BMGE527.0001	Balon fund rotund, 2 gaturi, 1000 ml	29/32	14/23	200	bucata
BMIE528.0001	Balon fund rotund, 2 gaturi, 1000 ml	29/32	19/26	200	bucata
BMFE529.0001	Balon fund rotund, 2 gaturi, 1000 ml	29/32	24/29	200	bucata
BMSE530.0001	Balon fund rotund, 2 gaturi, 2000 ml	24/29	19/26	231	bucata
BMVE531.0001	Balon fund rotund, 2 gaturi, 2000 ml	24/29	24/29	231	bucata
BMAE532.0001	Balon fund rotund, 2 gaturi, 2000 ml	29/32	19/26	234	bucata
BMCE533.0001	Balon fund rotund, 2 gaturi, 2000 ml	29/32	24/29	234	bucata

BALON CU FUND ROTUND, CU 3 GATURI

- Confectionat din sticla borosilicata 3,3 termorezistenta
- SN= diametrul partii de imbinare/ inaltimea partii de imbinare (mm)



Cod	Produs	SN-ul gatului central	SN-ul gaturilor laterale	H _{totala} (mm)	Mod de ambalare
BMRE534.0001	Balon fund rotund, 3 gaturi, 100 ml	19/26	14/23	106	bucata
BMBE535.0001	Balon fund rotund, 3 gaturi, 250 ml	19/26	14/23	134	bucata
BMDE536.0001	Balon fund rotund, 3 gaturi, 250 ml	19/26	19/26	134	bucata
BMZE537.0001	Balon fund rotund, 3 gaturi, 250 ml	24/29	14/23	139	bucata
BMPE538.0001	Balon fund rotund, 3 gaturi, 250 ml	24/29	19/26	139	bucata
BMLE539.0001	Balon fund rotund, 3 gaturi, 250 ml	24/29	24/29	139	bucata
BMBE540.0001	Balon fund rotund, 3 gaturi, 250 ml	29/32	14/23	145	bucata
BMGE541.0001	Balon fund rotund, 3 gaturi, 250 ml	29/32	19/26	145	bucata
BMXE542.0001	Balon fund rotund, 3 gaturi, 500 ml	19/26	14/23	158	bucata
BMQE543.0001	Balon fund rotund, 3 gaturi, 500 ml	19/26	19/26	158	bucata
BMTE544.0001	Balon fund rotund, 3 gaturi, 500 ml	24/29	14/23	161	bucata
BMNE545.0001	Balon fund rotund, 3 gaturi, 500 ml	24/29	19/26	161	bucata
BMSE546.0001	Balon fund rotund, 3 gaturi, 500 ml	24/29	24/29	161	bucata
BMBE547.0001	Balon fund rotund, 3 gaturi, 500 ml	29/32	14/23	168	bucata
BMPE548.0001	Balon fund rotund, 3 gaturi, 500 ml	29/32	19/26	167	bucata
BMLE549.0001	Balon fund rotund, 3 gaturi, 1000 ml	19/26	14/23	188	bucata
BMME550.0001	Balon fund rotund, 3 gaturi, 1000 ml	19/26	19/26	188	bucata
BMNE551.0001	Balon fund rotund, 3 gaturi, 1000 ml	24/29	14/23	196	bucata
BMCE552.0001	Balon fund rotund, 3 gaturi, 1000 ml	24/29	19/26	196	bucata
BMAE553.0001	Balon fund rotund, 3 gaturi, 1000 ml	24/29	24/29	196	bucata
BMDE554.0001	Balon fund rotund, 3 gaturi, 1000 ml	29/32	14/23	200	bucata
BMHE555.0001	Balon fund rotund, 3 gaturi, 1000 ml	29/32	19/26	200	bucata
BMKE556.0001	Balon fund rotund, 3 gaturi, 1000 ml	29/32	24/29	200	bucata
BMOE557.0001	Balon fund rotund, 3 gaturi, 1000 ml	29/32	29/32	200	bucata
BMWE558.0001	Balon fund rotund, 3 gaturi, 2000 ml	24/29	24/29	231	bucata
BMSE559.0001	Balon fund rotund, 3 gaturi, 2000 ml	29/32	19/26	234	bucata
BMBE560.0001	Balon fund rotund, 3 gaturi, 2000 ml	29/32	24/29	234	bucata
BMXE561.0001	Balon fund rotund, 3 gaturi, 2000 ml	29/32	29/32	234	bucata
BMPE562.0001	Balon fund rotund, 3 gaturi, 3000 ml	29/32	24/29	281	bucata

BMRE563.0001	Balon fund rotund, 3 gaturi, 5000 ml	40/38	24/29	333	bucata
BMTE564.0001	Balon fund rotund, 3 gaturi, 5000 ml	50/42	24/29	337	bucata

BALON CU FUND ROTUND, CU 4 GATURI

- Confectionat din sticla borosilicata 3,3 termorezistenta
- SN= diametrul partii de imbinare/ inaltimea partii de imbinare (mm)



Cod	Produs	SN-ul gatului central	SN-ul gaturilor laterale simetrice	SN-ul gatului lateral asimetric	H _{totala} (mm)	Mod de ambalare
BMYE565.0001	Balon fund rotund, 4 gaturi, 500 ml	24/29	19/26	14/23	161	bucata
BMFE566.0001	Balon fund rotund, 4 gaturi, 500 ml	24/29	19/26	19/26	161	bucata
BMZE567.0001	Balon fund rotund, 4 gaturi, 500 ml	24/29	19/26	14/23	161	bucata
BMCE568.0001	Balon fund rotund, 4 gaturi, 500 ml	24/29	24/29	14/23	161	bucata
BMLE569.0001	Balon fund rotund, 4 gaturi, 500 ml	29/32	19/26	14/23	167	bucata
BMUE570.0001	Balon fund rotund, 4 gaturi, 500 ml	29/32	24/29	14/23	167	bucata
BMAE571.0001	Balon fund rotund, 4 gaturi, 1000 ml	24/29	19/26	14/23	196	bucata

BMOE572.0001	Balon fund rotund, 4 gaturi, 1000 ml	24/29	19/26	19/26	196	bucata
BMFE573.0001	Balon fund rotund, 4 gaturi, 2000 ml	24/29	19/26	19/26	231	bucata
BMTE574.0001	Balon fund rotund, 4 gaturi, 2000 ml	29/32	19/26	14/23	234	bucata
BMRE575.0001	Balon fund rotund, 4 gaturi, 2000 ml	29/32	24/29	14/23	234	bucata
BMBE576.0001	Balon fund rotund, 4 gaturi, 3000 ml	29/32	24/29	14/23	281	bucata
BMEE577.0001	Balon fund rotund, 4 gaturi, 5000 ml	40/38	19/26	14/23	333	bucata
BMCE578.0001	Balon fund rotund, 4 gaturi, 5000 ml	40/38	24/29	14/23	333	bucata