

CAPSULE DIN PORTELAN

CAPSULA DE EVAPORARE

- Confectionata din portelan glazurat



Cod	Produs	Diametrul exterior	Inaltimea	Volum	Mod de ambalare
CRXC5147.0010	Capsula de evaporare	40 mm	16 mm	10 ml	10 bucati
CREC3028.0010	Capsula de evaporare	50 mm	20 mm	20 ml	10 bucati
CRGC3029.0010	Capsula de evaporare	63 mm	25 mm	40 ml	10 bucati
CRSC3030.0010	Capsula de evaporare	70 mm	26 mm	50 ml	10 bucati
CRZC3031.0010	Capsula de evaporare	80 mm	34 mm	75 ml	10 bucati
CRRC3032.0010	Capsula de evaporare	85 mm	36 mm	110 ml	10 bucati
CRTC3033.0001	Capsula de evaporare	100 mm	40 mm	180 ml	bucata
CRYC3034.0001	Capsula de evaporare	110 mm	45 mm	220 ml	bucata
CRIC3035.0001	Capsula de evaporare	115 mm	47 mm	310 ml	bucata
CRPC3036.0001	Capsula de evaporare	150 mm	45 mm	370 ml	bucata
CRBC3037.0001	Capsula de evaporare	160 mm	64 mm	640 ml	bucata
CRGC5148.0001	Capsula de evaporare	200 mm	80 mm	1000 ml	bucata

CRCC5149.0001	Capsula de evaporare	260 mm	86 mm	2200 ml	bucata
CRKC5150.0001	Capsula de evaporare	310 mm	110 mm	3400 ml	bucata