


## CRIC PENTRU LABORATOR

### CRIC DE LABORATOR

- Este confectionat din metal usor
- Foarte robust, poate sustine in pozitia initiala, pana la 60 de kg



Cod	Produs	Dimensiuni (L x l)	Inaltimea (ajustabila)	Greutate care poate fi ridicata	Mod de ambalare
CRDC3167.0001	Cric de laborator	110 x 100 mm	60 -200 mm	4 kg	bucata
CRFC3168.0001	Cric de laborator	160 x 130 mm	60 – 250 mm	7 kg	bucata
CRTC3169.0001	Cric de laborator	200 x 200 mm	70 – 260 mm	7 kg	bucata
CRGC3170.0001	Cric de laborator	200 x 230 mm	70 – 260 mm	7 kg	bucata
CRZC3171.0001	Cric de laborator	300 x 300 mm	130 – 465 mm	12 kg	bucata
CRXC3172.0001	Cric de laborator	400 x 400 mm	170 – 615 mm	15 kg	bucata

## CRIC DE LABORATOR

- Placa este confectionata din otel inoxidabil, grilajul mobil din aluminiu
- Prevazut cu picioare din cauciuc ce previn alunecarea
- Greutatea maxima suportata este de 7 kg



Cod	Produs	Dimensiuni (L x l)	Inaltime (ajustabila)	Greutate	Mod de ambalare
CRLC3173.0001	Cric de laborator	130 x 160 mm	60 -260 mm	1330 g	bucata
CRPC3174.0001	Cric de laborator	200 x 200 mm	60 – 265 mm	1970 g	bucata