

## EPRUBETE

### EPRUBETE FIOLEX

- In partea superioara prezinta un guler
- Rezistenta la schimbari bruste de temperatura si la incalzire locala



Cod	Produs	Diametrul exterior (mm)	Lungime (mm)	Grosimea peretelui (mm)	Mod de ambalare
ERIC379.0100	Eprubeta Fiolax 2 ml	8	70	0,4-0,5	100 bucati
ERFC380.0100	Eprubeta Fiolax 4 ml	10	75	0,4-0,5	100 bucati
ERO381.0100	Eprubeta Fiolax 6 ml	10	100	0,4-0,5	100 bucati
ERWC382.0100	Eprubeta Fiolax 6,5 ml	12	75	0,4-0,5	100 bucati
ERQC383.0100	Eprubeta Fiolax 9 ml	12	100	0,4-0,5	100 bucati
ERAC384.0100	Eprubeta Fiolax 16 ml	14	130	0,4-0,5	100 bucati
ERTC385.0100	Eprubeta Fiolax 20 ml	16	130	0,5-0,6	100 bucati
ERUC386.0100	Eprubeta Fiolax 25 ml	16	160	0,5-0,6	100 bucati
ERXC387.0100	Eprubeta Fiolax 35 ml	18	180	0,5-0,6	100 bucati
ERV388.0100	Eprubeta Fiolax 39 ml	20	150	0,5-0,6	100 bucati

ERN389.0100	Eprubeta Fiolax 45 ml	20	180	0,5-0,6	100 bucati
ERGC390.0050	Eprubeta Fiolax 60 ml	25	150	0,6-0,7	50 bucati
ERJC391.0050	Eprubeta Fiolax 80 ml	25	200	0,6-0,7	50 bucati
ERSC392.0050	Eprubeta Fiolax 110 ml	30	200	0,7-0,8	50 bucati

## EPRUBETE SIMPLE

- Sticla borosilicata de 3,3 mm termorezistenta
- Negradate



Cod	Produs	Øext (mm)	Inaltime (mm)	Grosime perete (mm)	Mod de ambalare
EMTE613.0001	Eprubeta simpla, 8/80 mm	8 ±0,5	80 ±1	0,8-1,0	bucata
EMPE614.0001	Eprubeta simpla, 10/100 mm	10 ±1	100 ±2	0,8-1,0	bucata
EMKE615.0001	Eprubeta simpla, 12/120 mm	12 ±1	120 ±2	0,8-1,0	bucata
EMLE616.0001	Eprubeta simpla, 14/140 mm	14 ±1	140 ±2	0,8-1,0	bucata
EMRE617.0001	Eprubeta simpla, 16/160 mm	16 ±1,2	160 ±3	0,8-1,0	bucata
EMSE618.0001	Eprubeta simpla, 18/180 mm	18 ±1,2	180 ±3	0,8-1,0	bucata

## EPRUBETE GRADATE

- Sticla borosilicata de 3,3 mm termorezistenta



Cod	Produs	Valoarea diviziunii (ml)	Dimensiuni ( $\varnothing_{\text{interior}} \times L$ ) mm	Mod de ambalare
EMOE309.0001	Eprubeta gradata, 5 ml	0,1	11 x 110	bucata
EMRE310.0001	Eprubeta gradata, 10 ml	0,2	14 x 130	bucata
EMVE311.0001	Eprubeta gradata, 15 ml	0,2	15 x 150	bucata
EMXE312.0001	Eprubeta gradata, 20 ml	0,5	17 x 165	bucata
EMWE313.0001	Eprubeta gradata, 25 ml	0,5	18 x 180	bucata
EMSE314.0001	Eprubeta gradata, 30 ml	0,5	19 x 180	bucata
EMJE315.0001	Eprubeta gradata, 50 ml	1,0	32 x 200	bucata

## EPRUBETE MINERALIZARE

- Sunt confectionate din sticla borosilicata 3,3 termorezistenta



Cod	Produs	Ø <sub>exterior</sub> x L	Grosime perete	Mod de ambalare
EMJE506.0001	Eprubeta mineralizare	42 x 300 mm	1,60 mm	bucata
EMDE507.0001	Eprubeta mineralizare	42 x 300 mm	3,20 mm	bucata